




acumed®

OsteoMed® ExtremiFix™ Cannulated Screw System

Mini & Small

Surgical Technique



Acumed® is a global leader of innovative orthopaedic and medical solutions.



We are dedicated to developing products, service methods, and approaches that improve patient care.



OsteoMed® ExtremiFix™ Cannulated Screw System | Mini & Small

The OsteoMed ExtremiFix Cannulated Screw System | Mini & Small is designed to provide a minimally invasive method of anatomical fixation in upper and lower extremities.

The ExtremiFix Cannulated Screw System | Mini & Small is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair, and fracture fixation of bones appropriate for the size of the device. Screws in the system are offered headed and headless in four diameters: 2.0, 2.4, 3.0, and 4.0 mm.

Refer to the provided Instructions for Use for the complete Indications, Contraindications, Warnings, and Instructions for Use including cleaning and sterilization details.

	Definition
Warning	Indicates critical information about a potential serious outcome to the patient or the user.
Caution	Indicates instructions that must be followed in order to ensure the proper use of the device.
Note	Indicates information requiring special attention.





Table of Contents

System Features	2
Instrument Overview	5
Surgical Technique Overview	6
Instrument and Screw Sizing	8
Surgical Techniques	9
Cannulated Screw Surgical Technique	9
Screw Removal Surgical Technique	13
Ordering Information	14

System Features


Cannulated Screws



System Features [continued]

Modular Tray Configurations

The ExtremiFix Cannulated Screw System | Mini & Small is designed as a modular storage system. Each module includes a complete set of instrumentation specific to that screw diameter. Outer cases contain general instrumentation and can store up to two screw modules. Additionally, independent ExtremiFix Cannulated Screw System | Mini & Small modules can supplement OsteoMed's Hand Plating System, ExtremiLock™ Foot Plating System, and ExtremiLock™ Ankle Plating System to treat a wider variety of upper and lower extremity applications.



System Features [continued]

Instrument Features



Double Sided Cannulated Depth Gauge
Dual-sided depth gauge provides screw length measurement for both 4" and 6" guide wires


Ratchet Handle
Quick release ratcheting handle with rubberized grip



Swivel Handle
Jewelers handle is designed for fine control and tactile feedback



Parallel Guide Wire
Allows surgeon to place two .045" guide wires perpendicular to the fracture



Ergonomic Design
Left and right-specific ratcheting medial malleolus clamp

Pointed Reduction Clamp
Designed to reduce and stabilize the medial malleolar fragment

Instrument Overview


General Instruments




Swivel Handle Driver, Quick Release
(323-1701)



Ratchet Handle, Quick Release
(323-1702)




Bone Clamp
(320-0102)



Bone Clamp, Curved Lobster Claw 5 1/2" Small
(316-0046)



Ø2.0/2.4/3.0/4.0 mm Screw Remover
(316-0039)



Ø2.0/2.4/3.0/4.0 mm Soft Tissue Protector/Wire Guide
(316-0008)



Ø2.0/2.4/3.0/4.0 mm Screw Extractor
(316-0009)



Cannulated Depth Gauge, Double Sided
(316-0036)

Optional Instruments*



Medial Malleolar Clamp, Left
(320-4045)



Medial Malleolar Clamp, Right
(320-4046)

*The Medial Malleolar Clamp (320-404x) is an optional instrument, and may not be required to complete installation of the ExtremiFix Cannulated Screw System | Mini & Small. For more information please see: 030-2231.

Instrument Overview [continued]

2.0/2.4 mm Mini Instruments



Ø1.7 mm Cannulated Drill, Short, Quick Release
(316-0513)



Ø1.7 mm Cannulated Drill, Long, Quick Release
(316-0514)



Ø2.0/2.4 mm Cannulated Countersink, Quick Release
(316-0511)



Ø2.0/2.4 mm Proximal Cortex Drill, Quick Release
(316-0516)



Ø1.7 mm Drill/Countersink Combo, Short, Quick Release
(316-0527)



Ø1.7 mm Drill/Countersink Combo, Long, Quick Release
(316-0528)



Ø2.0/2.4 mm T7 Driver, Cannulated, Short, Quick Release
(316-0544)




Ø2.0/2.4 mm T7 Driver, Cannulated, Long, Quick Release
(316-0545)


Surgical Technique Overview

Cannulated Screw Surgical Technique

Reduction and Provisional Fixation



Insert K-wire



Measure Screw Length



Instrument Overview [continued]

3.0/4.0 mm Small Instruments



Ø3.0/4.0 mm Cannulated Countersink, Quick Release
(316-0501)



Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release
(316-0554)



Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release
(316-0555)

3.0 mm Small Instruments



Ø2.3 mm Cannulated Drill, Quick Release
(316-0532)



Ø3.2 mm Proximal Cortex Drill, Quick Release
(316-0533)



Ø2.3 mm Drill/Countersink Combo, Short, Quick Release
(316-0537)



Ø2.3 mm Drill/Countersink Combo, Long, Quick Release
(316-0538)

4.0 mm Small Instruments



Ø2.7 mm Cannulated Drill, Quick Release
(316-0534)



Ø3.7 mm Proximal Cortex Drill, Quick Release
(316-0543)




Ø2.7 mm Drill/Countersink Combo, Short, Quick Release
(316-0547)



Ø2.7 mm Drill/Countersink Combo, Long, Quick Release
(316-0548)

Predrill (Recommended)



Countersink or Cortex Drill (Recommended)



Insert Screw



Screw Removal (Optional)



Instrument and Screw Sizing

Instrumentation and screw modules are color-coded based on screw diameter to aid in identification.



Diameter	2.0 mm Mini		2.4 mm Mini		3.0 mm Small		4.0 mm Small	
Material	Titanium Alloy							
Head Type	Headed	Headless	Headed	Headless	Headed	Headless	Headed	Headless
Length	6–42 mm	10–42 mm	6–50 mm	10–50 mm	10–40 mm	10–40 mm	12–52 mm	12–52 mm
K-wire	Ø0.035" x 4" and Ø0.035" x 6"				Ø0.045" x 4" and Ø0.045" x 6"			
Pilot Drill*	Ø1.7 mm				Ø2.3 mm		Ø2.7 mm	
Drive Type*	T7				T10			

*Association for Osteosynthesis (AO) standard quick release connection.

2.0 mm Mini Cannulated Screws



2.0 mm x 6–42 mm Cannulated Headed Screws
(318-20XX)
2 mm increments



2.0 mm x 10–42 mm Cannulated Headless Screws
(315-20XX)
2 mm increments



2.0 mm Washer
(316-0125)

2.4 mm Mini Cannulated Screws



2.4 mm x 6–50 mm Cannulated Headed Screws
(318-24XX)
2 mm increments



2.4 mm x 10–50 mm Cannulated Headless Screws
(315-24XX)
2 mm increments



2.4 mm Washer
(316-0126)

3.0 mm Small Cannulated Screws



3.0 mm x 10–40 mm Cannulated Headed Screws
(318-30XX)
2 mm increments




3.0 mm x 10–40 mm Cannulated Headless Screws
(315-30XX)
2 mm increments



3.0 mm Washer
(316-0225)

4.0 mm Small Cannulated Screws



4.0 mm x 12–52 mm Cannulated Headed Screws
(318-40XX)
2 mm increments up to 40 mm in Length
4 mm increments up to 52 mm in Length



4.0 mm x 12–52 mm Cannulated Headless Screws
(315-40XX)
2 mm increments up to 40 mm in Length
4 mm increments up to 52 mm in Length



4.0 mm Washer
(316-0226)

Foot & Ankle Applications

Headed and Headless Cannulated Screws

1 Reduction and Provisional Fixation

Using the preferred technique to reduce bone segments and provisionally fixate with a clamp or K-wire. A separate K-wire may provide increased rotational stability.




Figure 1

2 Insert K-wire

Place the Soft Tissue Protector/Wire Guide (316-0008) through a stab or formal open incision. Through the soft tissue protector, insert the K-wire (316-0XXX), across the fixation site to the appropriate depth (Figure 2).

Caution: Careful dissection and the soft tissue protector should be used during percutaneous placement.




Figure 2

K-wires for 2.0/2.4 mm Screws

316-0107	Ø.035 x 4" K-wire
316-0106	Ø.035 x 6" K-wire
316-0114	Ø.035 x 4" Partially Threaded K-wire
316-0119	Ø.035 x 4" Double Trocar K-wire
316-0117	Ø.035 x 6" Double Trocar K-wire

K-wires for 3.0/4.0 mm Screws

316-0003	Ø.045 x 4" K-wire
316-0004	Ø.045 x 6" K-wire
316-0006	Ø.045 x 6" Partially Threaded K-wire
316-0123	Ø.045 x 4" Double Trocar K-wire
316-0122	Ø.045 x 6" Double Trocar K-wire



Ø2.0/2.4/3.0/4.0 mm Soft Tissue Protector/Wire Guide (316-0008)



K-wire (316-0XXX)

Foot & Ankle Applications

Headed and Headless Cannulated Screws [continued]



Figure 3

3 Measure Screw Length

Slide the Cannulated Depth Gauge, Double Sided (316-0036) over the K-wire until the tip of the depth gauge contacts the surface of the bone (Figure 3). Make sure to use the side of the depth gauge that corresponds to the length of guide wire used.

Note: After screw length has been recorded and when appropriate, the guide wire may be advanced to maintain position during pre-drilling.




Figure 4

4 Predrill (Recommended)

The ExtremiFix Mini & Small cannulated screws are self-drilling and self-tapping.

Note: In cases of dense cortical bone, pre-drilling with the appropriately sized cannulated drill is recommended.

Drills for 2.0/2.4 mm Screws

316-0513	Ø1.7 mm Cannulated Drill, Short, Quick Release
316-0514	Ø1.7 mm Cannulated Drill, Long, Quick Release

Drill for 3.0/4.0 mm Screws

316-0532	Ø2.3 mm Cannulated Drill, Quick Release
----------	---

Use the selected Cannulated Drill (316-05XX) through the soft tissue protector to drill a pilot hole (Figure 4).

Caution: To maintain K-wire position, drilling should not go beyond the tip of the wire.



Cannulated Depth Gauge, Double Sided (316-0036)



K-wire (316-0XXX)



Ø1.7 mm Cannulated Drill, Short, Quick Release (316-0513)



Ø1.7 mm Cannulated Drill, Long, Quick Release (316-0514)



Ø2.3 mm Cannulated Drill, Quick Release (316-0532)

Foot & Ankle Applications

Headed and Headless Cannulated Screws [continued]

5 Countersink or Cortex Drill (Recommended)

Cortex Drills and Countersinks (316-05XX) should be placed over the K-wire to create a recess in the near cortex just deep enough to capture the screw head (Figure 5).

Caution: For Headed Screws, countersinking should be done by hand to avoid over-penetrating the proximal cortex.

Note: Countersinking may reduce stress risers and screw head prominence, particularly when guide wire placement is not perpendicular to the bone surface. For Headless Screws, Cortex Drills are available to create a pilot hole for proximal screw threads in dense cortical bone.




Figure 5

Countersink/Cortex Drill for 2.0/2.4 mm Screws

316-0511	Ø2.0/2.4 mm Cannulated Countersink, Quick Release
316-0516	Ø2.4 mm Proximal Cortex Drill

Countersink/Cortex Drill for 3.0/4.0 mm Screws

316-0501	Ø3.0/4.0 mm Cannulated Countersink, Quick Release
316-0533	Ø3.2 mm Proximal Cortex Drill, Quick Release



Ø2.0/2.4 mm
Cannulated Countersink,
Quick Release
(316-0511)



Ø2.4 mm Proximal
Cortex Drill
(316-0516)



Ø3.0/4.0 mm Cannulated
Countersink,
Quick Release
(316-0501)



Ø3.2 mm Proximal
Cortex Drill,
Quick Release
(316-0533)



K-wire
(316-0XXX)

Foot & Ankle Applications

Headed and Headless Cannulated Screws [continued]

Figure 6




Figure 7

6 Insert Screw

The Screw Remover (316-0039) may be helpful to remove the desired screw from the module. Verify screw length using the length gauge on the screw module (Figure 6).

Slide the selected Cannulated Screw (318-XXXX or 315-XXXX) over the K-wire.

Using the Cannulated Driver (316-05XXX) and Driver Handle (323-170X), drive the screw head until fully seated (Figure 7).

Driver for 2.0/2.4 mm Screws

316-0544	Ø2.0/2.4 mm T7 Driver, Cannulated, Short, Quick Release
316-0545	Ø2.0/2.4 mm T7 Driver, Cannulated, Long, Quick Release

Drivers for 3.0/4.0 mm Screws

316-0554	Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release
316-0555	Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release

It may be necessary to drive the screw further into the bone to prevent screw head prominence.







Verify final screw placement under fluoroscopy then remove and discard the K-wire.

Repeat the overall surgical technique, as necessary, for additional screw placement.

Washers

316-0125	Ø2.0 mm Washer
316-0126	Ø2.4 mm Washer
316-0225	Ø3.0 mm Washer
316-0226	Ø4.0 mm Washer

Note: To provide additional contact area with the bone, optional Washers (316-0XXX) may be placed over the K-wire, prior to a headed screw.



Screw Removal Technique [Optional]

1 Screw Removal

Locate the implant with intraoperative imaging.

Palpate the screw and remove surrounding soft tissue to gain maximum exposure.

Insert a K-wire (316-0XXX) through the screw cannula, if possible.

Engage screw with a Cannulated Driver (316-05XX) and rotate counterclockwise until the screw is removed (Figure 1).

If the screw cannot be removed with a cannulated driver, insert the Screw Extractor (316-0009) into the cannulation. The screw head does not need to be intact during extraction. As long as the screw canal is exposed, the tip of the screw extractor can be inserted into the cannula.

While applying axial force, rotate the screw extractor counterclockwise until the screw is removed.

Warning: It is important to keep the extractor in line with the screw to be removed. Failure to do so may cause the extractor to break.




Figure 1



Ø2.0/2.4/3.0/4.0 mm
Screw Extractor
(316-0009)



K-wire
(316-0XXX)



Swivel Handle Driver,
Quick Release
(323-1701)



Ratchet Handle,
Quick Release
(323-1702)



Ø2.0/2.4 mm T7 Driver,
Cannulated, Short,
Quick Release
(316-0544)



Ø2.0/2.4 mm T7 Driver,
Cannulated, Long,
Quick Release
(316-0545)



Ø3.0/4.0 mm T10
Driver, Cannulated,
Short, Quick Release
(316-0554)







Ø3.0/4.0 mm T10
Driver, Cannulated,
Long, Quick Release
(316-0555)

Ordering Information

Outer Trays					
Outer Trays - Mini System		Outer Trays - Small System			
1	Outer Tray Base, ExtremiFix Mini & Small	316-1501	3	Outer Tray Base, ExtremiFix Mini & Small	316-1501
2	Outer Tray Lid, ExtremiFix Mini & Small	316-1502	4	Outer Tray Lid, ExtremiFix Mini & Small	316-1502
Tray Components					
5	Inner Platter, Ratchet and Swivel Handles	316-1503	8	Inner Platter, Clamps and Common Instruments	316-1504
6	Swivel Handle Driver, Quick Release	323-1701	9	Bone Clamp	320-0102
7	Ratchet Handle, Quick Release	323-1702	10	Bone Clamp, Curved Lobster Claw 5½" Small	316-1501
			11	Ø2.0/2.4/3.0/4.0 mm Cannulated Screw Remover	316-0039
			12	Ø2.0/2.4/3.0/4.0 mm Soft Tissue Protector	316-0008
			13	Ø2.0/2.4/3.0/4.0 mm Screw Extractor	316-0009
Optional Instruments*					
14	Medial Malleolar Clamp, Left	320-4045	15	Medial Malleolar Clamp, Right	320-4046

*The Medial Malleolar Clamp is an optional instrument, and may not be required to complete installation of the ExtremiFix Cannulated Screw System | Mini & Small. For more information please see: 030-2231.

Note: To learn more about the full line of Acumed innovative surgical solutions, please contact your authorized Acumed distributor, call 888.627.9957, or visit www.acumed.net.



Ordering Information [continued]

Tray Components — 2.0 mm Screw Module

1	Ø2.0 mm Screw & Instrument Module, ExtremiFix Mini	316-1520	2	Ø1.7 mm Cannulated Drill, Short, Quick Release	316-0513
3	Ø1.7 mm Cannulated Drill, Long, Quick Release	316-0514	4	Ø2.0/2.4 mm Cannulated Countersink, Quick Release	316-0511
5	Ø2.0/2.4 mm Proximal Cortex Drill, Quick Release	316-0516	6	Ø1.7 mm Drill/Countersink Combo, Short, Quick Release	316-0527
7	Ø1.7 mm Drill/Countersink Combo, Long, Quick Release	316-0528	8	Ø2.0/2.4 mm T7 Driver, Cannulated, Short, Quick Release	316-0544
9	Ø2.0/2.4 mm T7 Driver, Cannulated, Long, Quick Release	316-0545			



Tray Components — 2.0 mm Screw Module

10	Ø2.0 mm Screw & Instrument Module, ExtremiFix Mini	316-1520	11	Cannulated Depth Gauge, Double sided 4"/6"	316-0036
12	Ø.035 x 6" K-wire	316-0106	13	Ø.035 x 6" Double Trocar K-wire	316-0117
14	Ø.035 x 4" K-wire	316-0107	15	Ø.035 x 4" Double Trocar K-wire	316-0119
16	Ø.035 x 4" Partially Threaded K-wire	316-0114	17	Ø2.0 mm Washer	316-0125

Tray Components — 2.0 mm Screw Module

18	2.0 mm Headed Cannulated Screws		19	2.0 mm Headless Cannulated Screws	
	2.0 x 6 mm Cannulated Headed Screw	318-2006		2.0 x 10 mm Cannulated Headless Screw	315-2010
	2.0 x 8 mm Cannulated Headed Screw	318-2008		2.0 x 12 mm Cannulated Headless Screw	315-2012
	2.0 x 10 mm Cannulated Headed Screw	318-2010		2.0 x 14 mm Cannulated Headless Screw	315-2014
	2.0 x 12 mm Cannulated Headed Screw	318-2012		2.0 x 16 mm Cannulated Headless Screw	315-2016
	2.0 x 14 mm Cannulated Headed Screw	318-2014		2.0 x 18 mm Cannulated Headless Screw	315-2018
	2.0 x 16 mm Cannulated Headed Screw	318-2016		2.0 x 20 mm Cannulated Headless Screw	315-2020
	2.0 x 18 mm Cannulated Headed Screw	318-2018		2.0 x 22 mm Cannulated Headless Screw	315-2022
	2.0 x 20 mm Cannulated Headed Screw	318-2020		2.0 x 24 mm Cannulated Headless Screw	315-2024
	2.0 x 22 mm Cannulated Headed Screw	318-2022		2.0 x 26 mm Cannulated Headless Screw	315-2026
	2.0 x 24 mm Cannulated Headed Screw	318-2024		2.0 x 28 mm Cannulated Headless Screw	315-2028
	2.0 x 26 mm Cannulated Headed Screw	318-2026		2.0 x 30 mm Cannulated Headless Screw	315-2030
	2.0 x 28 mm Cannulated Headed Screw	318-2028		2.0 x 32 mm Cannulated Headless Screw	315-2032
	2.0 x 30 mm Cannulated Headed Screw	318-2030		2.0 x 34 mm Cannulated Headless Screw	315-2034
	2.0 x 32 mm Cannulated Headed Screw	318-2032		2.0 x 36 mm Cannulated Headless Screw	315-2036
	2.0 x 34 mm Cannulated Headed Screw	318-2034		2.0 x 38 mm Cannulated Headless Screw	315-2038
	2.0 x 36 mm Cannulated Headed Screw	318-2036		2.0 x 40 mm Cannulated Headless Screw	315-2040
	2.0 x 38 mm Cannulated Headed Screw	318-2038		2.0 x 42 mm Cannulated Headless Screw	315-2042
	2.0 x 40 mm Cannulated Headed Screw	318-2040			
	2.0 x 42 mm Cannulated Headed Screw	318-2042			

Ordering Information [continued]



Ordering Information [continued]

Tray Components — 2.4 mm Screw Module

1	Ø2.4 mm Screw & Instrument Module, ExtremiFix Mini	316-1524	See previous spread for instrument callouts
----------	--	----------	--

Tray Components — 2.4 mm Screw Module

10	Ø2.4 mm Screw & Instrument Module, ExtremiFix Mini	316-1524	11	Cannulated Depth Gauge, Double sided 4"/6"	316-0036
12	Ø.035 x 6" K-wire	316-0106	13	Ø.035 x 6" Double Trocar K-wire	316-0117
14	Ø.035 x 4" K-wire	316-0107	15	Ø.035 x 4" Double Trocar K-wire	316-0119
16	Ø.035 x 4" Partially Threaded K-wire	316-0114	17	Ø2.4 mm Washer	316-0126

Tray Components — 2.4 mm Screw Module



18 2.4 mm Headed Cannulated Screws

2.4 x 6 mm Cannulated Headed Screw	318-2406
2.4 x 8 mm Cannulated Headed Screw	318-2408
2.4 x 10 mm Cannulated Headed Screw	318-2410
2.4 x 12 mm Cannulated Headed Screw	318-2412
2.4 x 14 mm Cannulated Headed Screw	318-2414
2.4 x 16 mm Cannulated Headed Screw	318-2416
2.4 x 18 mm Cannulated Headed Screw	318-2418
2.4 x 20 mm Cannulated Headed Screw	318-2420
2.4 x 22 mm Cannulated Headed Screw	318-2422
2.4 x 24 mm Cannulated Headed Screw	318-2424
2.4 x 26 mm Cannulated Headed Screw	318-2426
2.4 x 28 mm Cannulated Headed Screw	318-2428
2.4 x 30 mm Cannulated Headed Screw	318-2430
2.4 x 32 mm Cannulated Headed Screw	318-2432
2.4 x 34 mm Cannulated Headed Screw	318-2434
2.4 x 36 mm Cannulated Headed Screw	318-2436
2.4 x 38 mm Cannulated Headed Screw	318-2438
2.4 x 40 mm Cannulated Headed Screw	318-2440
2.4 x 42 mm Cannulated Headed Screw	318-2442
2.4 x 44 mm Cannulated Headed Screw	318-2444
2.4 x 46 mm Cannulated Headed Screw	318-2446
2.4 x 48 mm Cannulated Headed Screw	318-2448
2.4 x 50 mm Cannulated Headed Screw	318-2450

19 2.4 mm Headless Cannulated Screws

2.4 x 10 mm Cannulated Headless Screw	315-2410
2.4 x 12 mm Cannulated Headless Screw	315-2412
2.4 x 14 mm Cannulated Headless Screw	315-2414
2.4 x 16 mm Cannulated Headless Screw	315-2416
2.4 x 18 mm Cannulated Headless Screw	315-2418
2.4 x 20 mm Cannulated Headless Screw	315-2420
2.4 x 22 mm Cannulated Headless Screw	315-2422
2.4 x 24 mm Cannulated Headless Screw	315-2424
2.4 x 26 mm Cannulated Headless Screw	315-2426
2.4 x 28 mm Cannulated Headless Screw	315-2428
2.4 x 30 mm Cannulated Headless Screw	315-2430
2.4 x 32 mm Cannulated Headless Screw	315-2432
2.4 x 34 mm Cannulated Headless Screw	315-2434
2.4 x 36 mm Cannulated Headless Screw	315-2436
2.4 x 38 mm Cannulated Headless Screw	315-2438
2.4 x 40 mm Cannulated Headless Screw	315-2440
2.4 x 42 mm Cannulated Headless Screw	315-2442
2.4 x 44 mm Cannulated Headless Screw	315-2444
2.4 x 46 mm Cannulated Headless Screw	315-2446
2.4 x 48 mm Cannulated Headless Screw	315-2448
2.4 x 50 mm Cannulated Headless Screw	315-2450

Ordering Information [continued]



Ordering Information [continued]

Tray Components — 3.0 mm Driver and Drill Module

1	Ø3.0 mm Screw & Instrument Module, ExtremiFix Small	316-1530	2	Ø2.3 mm Cannulated Drill, Quick Release	316-0532
3	Ø3.0/4.0 mm Cannulated Countersink, Quick Release	316-0501	4	Ø3.2 mm Proximal Cortex Drill, Quick Release	316-0533
5	Ø2.3 mm Drill/Countersink Combo, Short, Quick Release	316-0537	6	Ø2.3 mm Drill/Countersink Combo, Long, Quick Release	316-0538
7	Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release	316-0554	8	Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release	316-0555



Tray Components — 3.0 mm Screw Module

9	Ø3.0 mm Screw & Instrument Module, ExtremiFix Small	316-1530	10	Cannulated Depth Gauge, Double sided 4"/6"	316-0036
11	Ø.045 x 6" K-wire	316-0004	12	Ø.045 x 6" Partially Threaded K-wire	316-0006
13	Ø.045 x 6" Double Trocar K-wire	316-0122	14	Ø.045 x 4" K-wire	316-0003
15	Ø.045 x 4" Double Trocar K-wire	316-0123	16	Ø3.0 mm Washer	316-0225

Tray Components — 3.0 mm Screw Module

17	3.0 mm Headed Cannulated Screws	18	3.0 mm Headless Cannulated Screws	
	3.0 x 10 mm Cannulated Headed Screw	318-3010	3.0 x 10 mm Cannulated Headless Screw	315-3010
	3.0 x 12 mm Cannulated Headed Screw	318-3012	3.0 x 12 mm Cannulated Headless Screw	315-3012
	3.0 x 14 mm Cannulated Headed Screw	318-3014	3.0 x 14 mm Cannulated Headless Screw	315-3014
	3.0 x 16 mm Cannulated Headed Screw	318-3016	3.0 x 16 mm Cannulated Headless Screw	315-3016
	3.0 x 18 mm Cannulated Headed Screw	318-3018	3.0 x 18 mm Cannulated Headless Screw	315-3018
	3.0 x 20 mm Cannulated Headed Screw	318-3020	3.0 x 20 mm Cannulated Headless Screw	315-3020
	3.0 x 22 mm Cannulated Headed Screw	318-3022	3.0 x 22 mm Cannulated Headless Screw	315-3022
	3.0 x 24 mm Cannulated Headed Screw	318-3024	3.0 x 24 mm Cannulated Headless Screw	315-3024
	3.0 x 26 mm Cannulated Headed Screw	318-3026	3.0 x 26 mm Cannulated Headless Screw	315-3026
	3.0 x 28 mm Cannulated Headed Screw	318-3028	3.0 x 28 mm Cannulated Headless Screw	315-3028
	3.0 x 30 mm Cannulated Headed Screw	318-3030	3.0 x 30 mm Cannulated Headless Screw	315-3030
	3.0 x 32 mm Cannulated Headed Screw	318-3032	3.0 x 32 mm Cannulated Headless Screw	315-3032
	3.0 x 34 mm Cannulated Headed Screw	318-3034	3.0 x 34 mm Cannulated Headless Screw	315-3034
	3.0 x 36 mm Cannulated Headed Screw	318-3036	3.0 x 36 mm Cannulated Headless Screw	315-3036
	3.0 x 38 mm Cannulated Headed Screw	318-3038	3.0 x 38 mm Cannulated Headless Screw	315-3038
	3.0 x 40 mm Cannulated Headed Screw	318-3040	3.0 x 40 mm Cannulated Headless Screw	315-3040

Ordering Information [continued]



Ordering Information [continued]

Tray Components — 4.0 mm Driver and Drill Module

1	Ø4.0 mm Screw & Instrument Module, ExtremiFix Small	316-1540	2	Ø2.7 mm Cannulated Drill, Quick Release	316-0534
3	Ø3.0/4.0 mm Cannulated Countersink, Quick Release	316-0501	4	Ø3.7 mm Proximal Cortex Drill, Quick Release	316-0543
5	Ø2.7 mm Drill/Countersink Combo, Short, Quick Release	316-0547	6	Ø2.7 mm Drill/Countersink Combo, Long, Quick Release	316-0548
7	Ø3.0/4.0 mm T10 Driver, Cannulated, Short, Quick Release	316-0554	8	Ø3.0/4.0 mm T10 Driver, Cannulated, Long, Quick Release	316-0555



Tray Components — 4.0 mm Screw Module

9	Ø4.0 mm Screw & Instrument Module, ExtremiFix Small	316-1540	10	Cannulated Depth Gauge, Double sided 4"/6"	316-0036
11	Ø.045 x 6" K-wire	316-0004	12	Ø.045 x 6" Partially Threaded K-wire	316-0006
13	Ø.045 x 6" Double Trocar K-wire	316-0122	14	Ø.045 x 4" K-wire	316-0003
15	Ø.045 x 4" Double Trocar K-wire	316-0123	16	Ø4.0 mm Washer	316-0226

Tray Components — 4.0 mm Screw Module

17	4.0 mm Headed Cannulated Screws	18	4.0 mm Headless Cannulated Screws	
	4.0 x 12 mm Cannulated Headed Screw	318-4012	4.0 x 12 mm Cannulated Headless Screw	315-4012
	4.0 x 14 mm Cannulated Headed Screw	318-4014	4.0 x 14 mm Cannulated Headless Screw	315-4014
	4.0 x 16 mm Cannulated Headed Screw	318-4016	4.0 x 16 mm Cannulated Headless Screw	315-4016
	4.0 x 18 mm Cannulated Headed Screw	318-4018	4.0 x 18 mm Cannulated Headless Screw	315-4018
	4.0 x 20 mm Cannulated Headed Screw	318-4020	4.0 x 20 mm Cannulated Headless Screw	315-4020
	4.0 x 22 mm Cannulated Headed Screw	318-4022	4.0 x 22 mm Cannulated Headless Screw	315-4022
	4.0 x 24 mm Cannulated Headed Screw	318-4024	4.0 x 24 mm Cannulated Headless Screw	315-4024
	4.0 x 26 mm Cannulated Headed Screw	318-4026	4.0 x 26 mm Cannulated Headless Screw	315-4026
	4.0 x 28 mm Cannulated Headed Screw	318-4028	4.0 x 28 mm Cannulated Headless Screw	315-4028
	4.0 x 30 mm Cannulated Headed Screw	318-4030	4.0 x 30 mm Cannulated Headless Screw	315-4030
	4.0 x 32 mm Cannulated Headed Screw	318-4032	4.0 x 32 mm Cannulated Headless Screw	315-4032
	4.0 x 34 mm Cannulated Headed Screw	318-4034	4.0 x 34 mm Cannulated Headless Screw	315-4034
	4.0 x 36 mm Cannulated Headed Screw	318-4036	4.0 x 36 mm Cannulated Headless Screw	315-4036
	4.0 x 38 mm Cannulated Headed Screw	318-4038	4.0 x 38 mm Cannulated Headless Screw	315-4038
	4.0 x 40 mm Cannulated Headed Screw	318-4040	4.0 x 40 mm Cannulated Headless Screw	315-4040
	4.0 x 44 mm Cannulated Headed Screw	318-4044	4.0 x 44 mm Cannulated Headless Screw	315-4044
	4.0 x 48 mm Cannulated Headed Screw	318-4048	4.0 x 48 mm Cannulated Headless Screw	315-4048
	4.0 x 52 mm Cannulated Headed Screw	318-4052	4.0 x 52 mm Cannulated Headless Screw	315-4052

Ordering Information [continued]





www.acumed.net

Acumed USA Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
Office: +1.888.627.9957
Office: +1.503.627.9957

OsteoMed USA Campus
3885 Arapaho Road
Addison, TX 75001
Office: +1.800.456.7779
Worldwide: 001.972.677.4600

Acumed Iberica Campus
C. de Álvaro Caballero, 14,
28023 Madrid, Spain
Office: +34.913.51.63.57

030-2218 Rev.B | Effective: 2022/03 | © 2022 OsteoMed®

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

OsteoMed LLC is a wholly owned subsidiary of Acumed LLC.

OsteoMed®, ExtremiLock™, and ExtremiFix™ are registered trademarks of OsteoMed LLC.
Acumed® is registered trademark of Acumed LLC.