## Biotrak® Resorbable Fixation System

Acumed® LLC  
5885 NE Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1.503.627.9957  
acumed.net

OBSOLETE  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for latest version.

77013001-N  
EFFECTIVE 01-2020

## TABLE OF CONTENTS

Click on a language to navigate to the page

English – US .....	3
Español – ES .....	45
Português – PT .....	89

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for latest version.

## ACUMED® BIOTRAK® RESORBABLE FIXATION SYSTEM

### FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON AND SUPPORTING HEALTHCARE PROFESSIONALS

**DESCRIPTION:** The ACUMED resorbable fixation devices are designed to provide fixation of specific fractures, osteotomies, and arthrodeses while they heal.

**INDICATIONS:** This fixation device is indicated for small bone fractures, osteotomies, and arthrodeses, cancellous fragments, and osteochondral fragments in the upper and lower extremities.

**CONTRAINDICATIONS:** This device is contraindicated in the presence of active or latent infection, sepsis, osteoporosis, insufficient quantity and/or quality of bone. Patients who are unwell or incapable of following postoperative care instructions are contraindicated for these devices. This device is not intended for use in any indication not covered above.

**IMPLANT MATERIAL SPECIFICATIONS:** The implants are made of PLLA (Poly-L-lactic acid).

**SURGICAL INSTRUMENT MATERIAL SPECIFICATIONS:** The instruments are made of various grades of stainless steel, titanium, aluminum, and polymers evaluated for biocompatibility.

**IMPLANT INFORMATION FOR USE:** The surgeon must select the type and size of implant that best meets the patient's surgical needs. Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

**SURGICAL INSTRUMENT INFORMATION FOR USE:** Instruments provided with this system may be single use or reusable.

- The user must refer to the instrument's label to determine whether the instrument is single use or reusable. Single use instruments are labeled with a "do not re-use" symbol as described in the Symbol Legend section, below.
- Single use instruments must be discarded after a single use.
- Reusable instruments have a limited lifespan. Prior to and after each use, reusable instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to drivers, drill bits and instruments used for cutting or implant insertion.

**SURGICAL TECHNIQUES:** Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be

familiar with relevant publications and consult with experienced associates regarding the procedure before use. Surgical Techniques can be found on the Acumed website (acumed.net).

**IMPLANT WARNINGS:** For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant, material, methods of application, instruments, and recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, and/or excessive activity. Device breakage and/or damage can occur when the implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device failing due to aforementioned causes. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. In addition, because of the thermal sensitivity of bioabsorbable materials, the device shall not be utilized if the dot in the middle of the temperature sticker has turned black. A black dot on the sensor signifies that the environmental temperature may have exceeded the softening temperature of the bioabsorbable material during storage and/or transit. Exceeding the softening temperature of the material can lead to degradation of the mechanical properties and/or warping.

**SURGICAL INSTRUMENT WARNINGS:** For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application and the recommended surgical technique. Instrument breakage or damage, as well as tissue damage, can occur when an instrument is subjected to excessive loads, excessive speeds, dense bone, improper use or unintended use. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments.

**IMPLANT PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect implants against scratching and nicking. Such stress concentrations can lead to failure. Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons. The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

**SURGICAL INSTRUMENT PRECAUTIONS:** Single use surgical instruments shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect instruments against scratching and nicking, as such stress concentrations can lead to failure.

**ADVERSE EFFECTS:** Fracture of the implant due to excessive loading, incomplete or inadequate healing, or excessive force during insertion; Implant migration and/or loosening; Sterile

inflammation as a result of a body reaction to the degradation products of the absorbable material; Pain, discomfort, or abnormal sensations due to the presence of an implant; Nerve damage resulting from surgical trauma; Bone necrosis or bone resorption.

## CLEANING:

**Implant Cleaning:** Implants are provided sterile only and should not be re-cleaned.

**Instrument Cleaning:** Acumed Instruments and Accessories must be thoroughly cleaned before reuse, following the guidelines below :

### Warnings & Precautions

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/ reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.

- Mineral oil or silicone lubricants should not be used on Acumed instruments.
- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralised and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

### Manual Cleaning/Disinfection Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood

the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.

3. Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
6. Inspect instruments under normal lighting for the removal of visible soil.
7. If visible soil is seen, repeat the sonication and rinse steps above.
8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

OBsolete  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for latest version.

### Combination Manual/Automated Cleaning and Disinfecting Instructions

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.
2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. *Note: Use of a sonicator will aid in thorough cleaning of instruments. Using a syringe or water jet will improve flushing of difficult to reach areas and any closely mated surface.*
3. Remove instruments from enzyme solution and rinse in deionized water for a minimum of one (1) minute.
4. Place instruments in a suitable washer/ disinfecter basket and process through a standard washer/disinfecter cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

Step	Description
1	Two (2) minute prewash with cold tap water
2	Twenty (20) second enzyme spray with hot tap water
3	One (1) minute enzyme soak
4	Fifteen (15) second cold tap water rinse (X2)
5	Two (2) minute detergent wash with hot tap water (146-150°F / 64-66°C)
6	Fifteen (15) second hot tap water rinse
7	Ten (10) second purified water rinse with optional lubricant (146-150°F / 64-66°C)
8	Seven (7) minute hot air dry (240°F / 116°C)
<i>Note: Follow washer/disinfecter manufacturer's instructions explicitly.</i>	

### Automated Cleaning/Disinfection Instructions

- Automated washer/dryer systems are not recommended as the only cleaning method for surgical instruments.
- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

## STERILITY:

System components may be provided sterile or nonsterile.

**Sterile Product:** Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Acumed does not recommend resterilization of sterile-packaged product. If sterile packaging is damaged, the incident must be reported to Acumed. The product must not be used, and must be returned to Acumed.

**Non-Sterile Product:** Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all ~~implants and~~ instruments must be considered nonsterile, and sterilized by the hospital prior to use.

- Follow current AORN “Recommended Practices for Sterilization in Perioperative Practice Settings” and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Consult your equipment manufacturer’s written instructions for specific sterilizer and load configuration instructions.
- **Gravity-displacement steam sterilization is not recommended.**
- Flash sterilization is not recommended.

### Steam Sterilizer Parameters<sup>1</sup>

<b>Storage case part numbers:</b>		Tray: 80-0485, 80-0517 Case lid: 80-0486, 80-0518
Gravity Displacement Steam Sterilization:		
NOT RECOMMENDED		
Pre-Vacuum Steam Sterilization:		
Condition:	Wrapped*	
Exposure Temperature:	270°F (132°C)	
Exposure Time:	4 minutes	
Dry Time:	30 minutes	
Cooling Time <sup>2</sup>	30 minutes	

<sup>1</sup>Values in this table reflect the minimum parameters validated to achieve the required Sterility Assurance Level (SAL), for a fully loaded tray with all parts placed appropriately.

<sup>2</sup>Cooling time describes the validated interval following dry time and prior to handling. This interval is included for safe handling and the prevention of contamination; see ANSI/AAMI ST79:2010 Section 8.8.1.

\* US customers shall use an FDA-cleared sterilization wrap.

**STORAGE INSTRUCTIONS:** Store implants in a cool dry place and keep away from direct sunlight. The implant should be stored in a temperature and humidity controlled environment below 120°F. Prior to use, inspect packaging for signs of tampering or water contamination. Check the expiration date on the box to ensure that the shelf life of the product has not been exceeded. Also, please inspect the temperature sensor on the box to ensure that the center has not turned black. If the temperature of the packaging has exceeded the softening temperature of the material during its lifetime, the central dot on the temperature sensor will be black and this device must not be implanted. Use the oldest lots first because the bioabsorbable material has a finite shelf life.

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Catalogue number
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. For Professional Use Only.**

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