ACUMED®

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STERILITY: This product is provided sterile and was exposed to a minimum dose of 25.0 kGy gamma irradiation. With the exception of the Polarus Cap, resterilization of the other rods and screws of this product line may only be performed if the original sterile package has been opened in error using one of the following methods.

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- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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POLARUS HUMERAL RODS, ACUMED ULNA, RADIUS, and FIBULA RODS

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The Acumed intramedullary rods and screws are designed to provide fixation of humeral, forearm, and fibula fractures while they heal.

INFORMATION FOR USE: Physiological dimensions limit the sizes of implant appliances. The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

INDICATIONS: Polarus Rods address fracture fixation and osteotomies of the humerus. Acumed Radius and Ulna Rods address fracture fixation and osteotomies of the forearm. Acumed Fibula Rods address fracture fixation of the fibula.

CONTRAINDICATIONS: Active or latent infection. Osteoporosis, insufficient quantity or quality of bone/soft tissue. Material sensitivity. If suspected, tests are performed prior to implantation. Sepsis. Patients who are unwilling or incapable of following postoperative care instructions. This device is not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the

possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility the device failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

ADVERSE EFFECTS: Fracture of the implant due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material. Pain, discomfort, or abnormal sensations due to the presence of an implant. Nerve damage resulting from the surgical trauma. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing.

STERILITY: This product is provided sterile and was exposed to a minimum dose of 25,000-Rads gamma irradiation. With the exception of the Polarus Cap Screw, resterilization of the other rods and screws of this product line may only be performed if the original sterile package has been opened in error using one of the following methods.

Sterilization Methods

Gravity Displacement Autoclave	Full Cycle Time: 270° F (132° C) for 30 minutes Dry Time: 65 minutes
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 14 minutes Dry Time: 65 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.

- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.