# ACUMED®

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## ACUMED REUSABLE SURGICAL INSTRUMENTS AND ACCESSORIES

### FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** Acumed Reusable Surgical Instruments and Accessories are available for a wide variety of surgical techniques and implant systems. Proper care and inspection after each use is critical.

**INFORMATION FOR USE:** Many of Acumed Instruments and Accessories are reusable; however they have a limited life span. Prior to and after each use, the instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to hex drivers, drill bits and instruments used for cutting or implant insertion.

**WARNINGS:** For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage can occur when an instrument is subjected to excessive loads, speeds, or dense bone. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments.

**CLEANING:** Acumed Instruments and Accessories must be thoroughly cleaned before reuse. Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Excess blood or debris should be wiped off to prevent it from drying onto the surface. Using an enzymatic cleaning product such as Enzol, prepare the cleaning solution according to manufacturer's instructions.

**NOTE:** For exposed springs, coils, or flexible features - Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface

with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.

**NOTE:** Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

Avoiding prolonged soaking, submerge the device in the enzymatic cleaning solution for 5 minutes. After 5 minutes, use a soft bristle brush to remove any detergent residue. Instruments with lumens also should be held under water (large side up) ensuring the lumen is flushed with water. Place the device in an ultrasonic cleaner filled with deionized water, ensuring that there is no air trapped in any of the crevices by gently agitating the device. Sonicate the device for 5 minutes. Remove the device from the ultrasonic cleaner and rinse in water by submerging it completely. Gently wipe the device with a clean, dry cloth, then air dry.

**NOTE:** Even surgical instruments manufactured from high grade stainless steel must be dried thoroughly to prevent rust formation, and all devices must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.

**STERILITY:** This product is provided non-sterile. Sterilization may be performed using one of the following methods.

#### Sterilization Methods

Gravity Displacement Autoclave	<b>Full Cycle Time:</b> 270° F (132° C) for 60 minutes <b>Dry Time:</b> 125 minutes
Pre-Vacuum Autoclave	<b>Full Cycle Time:</b> 270° F (132° C) for 25 minutes <b>Dry Time:</b> 70 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

**STORAGE INSTRUCTIONS:** Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination.

**CAUTION:** For Professional Use Only.

**CAUTION:** Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.

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**WARNINGS:** For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage can occur when an instrument is subjected to excessive loads, speeds, or dense bone. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments

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Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.

- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2006 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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