ACUMED®

BONE PLATE SYSTEMS

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed® LLC
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-55-C
Effective 10/2013

STERILE R

REF

LOT

CE

0473

* P K G I - 5 5 - C *

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

KNOGLEPLADESYSTEMER

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healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment.

1. Shellock, F., ed. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking, as such stress concentrations can lead to failure.

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Sterilization Methods

Gravity Displacement Autoclave	NOT RECOMMENDED
Pre-Vacuum Autoclave	Full Cycle Time: 132° C (270° F) for 4 minutes Dry Time: 30 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN® "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest first.

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¹ Shellock, F. G. *Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition*. Biomedical Research Publishing Group, 2011.

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¹ Shellock, J.G., *Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition*, Biomedical Research Publishing Group, 2011.

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CONGRUENT BONE PLATE SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The Acumed Congruent Bone Plate System, Acumed Lower Extremity Bone Plate System, Polarus® Locking Proximal Humeral Plate System, Scapula Congruent Bone Plate System of bone plates, screws and accessories are designed to provide fixation for fractures, fusions or osteotomies.

INFORMATION FOR USE: Physiological dimensions limit the sizes of implant appliances. The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

INDICATIONS: The Acumed Congruent Bone Plate System includes plates, screws and accessories designed to provide fixation during fractures, fusions, and osteotomies for the clavicle, humerus, radius, ulna, metacarpal, metatarsal, malleolus, tibia and fibula. The Acumed Lower Extremity Bone Plate System includes plates, screws and accessories designed to provide fixation during fractures, fusions, and osteotomies for the phalanges, metatarsals, tarsals fibula, and tibia bones. The Acumed Polarus® Locking Proximal Humeral System includes plates, screws and accessories designed for the Humerus. The Acumed Scapula Congruent Bone Plate System provides fixation during fractures, fusions, and osteotomies for the scapula.

CONTRAINDICATIONS: Contraindications for the system are active or latent infection; sepsis; osteoporosis, insufficient quantity or quality of bone/soft tissue; and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete

healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment¹.

¹ Sherlock, K.G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Precious stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking, as such stress concentrations can lead to failure.

ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

STERILITY: The product is provided sterile and non-sterile. Sterile device was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a sterile device may only be performed if the original sterile package has been opened in error. Sterilization may be performed using one of the following methods.

Sterilization Methods

Sterilization Method	Recommendation
Gravity Displacement Autoclave	NOT RECOMMENDED
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 4 minutes Dry Time: 30 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN® Recommended Practices for Sterilization in Perioperative Practice Settings® and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

CAUTION: For Professional Use Only.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.