 **acumed**[®]
Multi-System Implant Components

CE
2797

MediMark[®] Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed[®] LLC
5885 NE Cornelius Pass Road
Hillsboro, OR 97124-9432
+1.503.627.9957
acumed.net

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

PKGI-73-H
EFFECTIVE 03-2019

TABLE OF CONTENTS

Click on a language to navigate to the page

| | |
|-----------------------|----|
| English – US | 3 |
| Dansk – DA | 8 |
| English – EN | 13 |
| Deutsch – DE | 18 |
| Ελληνικά – EL | 24 |
| Español – ES | 30 |
| Français – FR | 36 |
| Italiano – IT | 42 |
| Nederlands – NL | 48 |
| Norsk – NO | 54 |
| Português – PT | 59 |
| Suomi – FI | 65 |
| Svenska – SV | 70 |
| Türkçe – TR | 75 |

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

ACUMED® MULTI-SYSTEM IMPLANT COMPONENTS

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: Acumed Multi-System Implant Components are implants utilized across multiple Acumed solution systems.

INDICATIONS: See the specific system's instructions for use to determine the indications for multi-system implant components within the context of a specific system. Specific system IFUs are assigned to the unique implant or tray components of a system.

CONTRAINDICATIONS: See the specific system's instructions for use to determine the contraindications for multi-system implant components within the context of a specific system. Specific system IFUs are assigned to the unique implant or tray components of a system.

IMPLANT MATERIAL SPECIFICATIONS: The implants are made of titanium alloy per ASTM F136 or cobalt-chrome per ASTM F1537.

IMPLANT INFORMATION FOR USE: The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

SURGICAL TECHNIQUES: Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use. Surgical Techniques can be found on the Acumed website (acumed.net).

IMPLANT WARNING S: For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/ or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to

delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/ or treatment to fail. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment ¹.

¹ Shellock, F. G. *Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition*. Biomedical Research Publishing Group, 2011.

IMPLANT PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Protect implants against scratching and nicking. Such stress concentrations can lead to failure. Bending plates multiple times may weaken the device and could lead to premature implant fracture and failure. Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons. The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue

damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

CLEANING:

Implant Cleaning: Implants should not be reused. Acumed does not recommend cleaning of implants provided sterile. Implants provided non-sterile that have not been used, but have become soiled, should be processed according to the following:

Warnings & Precautions

- Sterilization of the implants should not be performed if the implant comes into contact with contamination (e.g. biological tissue contact, such as bodily fluids/ blood) unless the single use device (SUD) has been reprocessed by an authorized facility who has received appropriate regulatory clearance for such. Cleaning a SUD after it comes into contact with human blood or tissue constitutes reprocessing.
- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- Users should wear appropriate personal protective equipment (PPE).

- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines, standards and hospital policies.

Manual Processing

Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

- Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
- Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
- Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
- Dry the implant using a clean, soft, lint-free cloth to avoid scratching the surface.

Ultrasonic Processing

Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.

- Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
- Clean implants ultrasonically for a minimum of 15 minutes.

- Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
- Dry the implant using a clean, soft, lint-free cloth to avoid scratching the surface.

Mechanical Processing

Equipment: Washer/disinfector, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

| Cycle | Minimum Time (minutes) | Minimum Temperature/Water | Type of Detergent |
|-------------|------------------------|--------------------------------------|-------------------------------|
| Pre-wash | 2 | Cold tap water | N/A |
| Enzyme Wash | 2 | Warm tap water | Neutral enzymatic pH ≤ 8.5 |
| Wash II | 5 | Warm tap water (>40°C) | Detergent with pH ≤ 8.5 |
| Rinse | 2 | Warm DI or purified water (>40°C) | N/A |
| Dry | 40 | 90°C | N/A |

STERILITY:

System components may be provided sterile or nonsterile.

Sterile Product: Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Acumed does not recommend resterilization of sterile-packaged product. If sterile packaging is damaged, the incident must be reported to Acumed. The product must not be used, and must be returned to Acumed.

Non-Sterile Product: Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Sterilization of nonsterile devices has been validated using the sterilization parameters listed below, where devices are provided in fully-loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79:

2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities

| | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| See a specific system's Instruction For Use to determine if shorter sterilization times are applicable for that system. Specific system IFUs are assigned to the unique implant or tray components of a system. | |
| Gravity Displacement Autoclave: NOT RECOMMENDED | |
| Pre-Vacuum Autoclave: | |
| Condition: | Wrapped |
| Exposure Temperature: | 270° F (132° C) |
| Exposure Time: | 4 minutes |
| Dry Time: | 30 minutes |

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

APPLICABILITY: These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

FURTHER INFORMATION: To request further material, please see the contact information listed on this document.

| SYMBOL LEGEND | |
|------------------------------------------------------------------------------------|-----------------------------------------------------|
|  | Consult instructions for use |
|  | Caution |
|  | Sterilized using ethylene oxide |
|  | Sterilized using irradiation |
|  | Use-by date |
|  | Catalogue number |
|  | Batch code |
|  | Authorized representative in the European Community |
|  | Manufacturer |
|  | Date of manufacture |
|  | Do not re-sterilize |
|  | Do not re-use |
|  | Upper limit of temperature |

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. For Professional Use Only.

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INDICATIONS: See the specific system's instructions for use to determine the indications for multi-system implant components within the context of a specific system. Specific system IFUs are assigned to the unique implant or tray components of a system .

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¹ *Shellock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.*

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1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
2. Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean, soft, lint-free cloth to avoid scratching the surface.

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1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
2. Clean implants ultrasonically for a minimum of 15 minutes.

3. Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
4. Dry the implant using a clean, soft, lint-free cloth to avoid scratching the surface.

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| Cycle | Minimum Time (minutes) | Minimum Temperature/Water | Type of Detergent |
|------------|------------------------|--------------------------------------|-------------------------------|
| Pre-wash | 2 | Cold tap water | N/A |
| EnzymeWash | 2 | Warm tap water | Neutral enzymatic pH ≤ 8.5 |
| Wash II | 5 | Warm tap water (>40°C) | Detergent with pH ≤ 8.5 |
| Rinse | 2 | Warm DI or purified water (>40°C) | N/A |
| Dry | 40 | 90°C | N/A |

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See a specific system's Instructions For Use to determine if shorter sterilization times are applicable for that system. Specific system IFUs are assigned to the unique implant or tray components of a system.

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