Sterilization In Aesculap® Rigid Sterilization Containers

US	Supplemental Information	2
DA	Supplerende Oplysninger.....	8
DE	Ergänzende Informationen.....	14
EL	Συμπληρωματικές Πληροφορίες	20
EN	Supplemental Information	26
ES	Información Adicional.....	32
FI	Täydentävät Tiedot	38
FR	Informations Complémentaires	44
IT	Informazioni Supplementari.....	50
NB	Tilleggsinformasjon	56
NL	Aanvullende Informatie	62
PT	Informação Adicional.....	68
SV	Kompletterande Information.....	74
TR	EK Bilgi	80

Visit www.acumed.net/ifu for the latest version.

Supplemental Information

Sterilization In Aesculap® Rigid Sterilization Containers

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The US instructions are intended for users in the United States and its territories.

Rx only

DESCRIPTION

The Acumed system configurations listed below have been validated for sterilization in the specified Aesculap sterilization containers.

INFORMATION FOR USE

Acumed Instructions for Use for implant and instrument systems can be found by part number at <http://www.acumed.net/ifu>. These instructions include product descriptions, information for use, indications, contraindications, warnings, precautions, adverse effects, cleaning instructions, and storage instructions for the personal attention of the operating surgeon.

Additional Aesculap Instructions for Use are provided by Aesculap. These instructions can be found by Aesculap Item Number at <http://www.aesculapusa.com> and must be reviewed prior to use of Aesculap sterile containers. Item Numbers can be found in the table below.

STERILITY

- Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately. Sterilization must be completed in adherence to the additional Instructions for Use provided by Aesculap.

CLEANING

- Cleaning of Aesculap sterilization containers should be completed in adherence to the Instructions for Use provided by Aesculap.
- Cleaning of the Acumed products should be completed in adherence to the part or system Instructions for Use provided by Acumed by product part number.

STERILIZATION

- These products are non-sterile, and must be sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately. Sterilization must be completed in adherence to the additional Instructions for Use provided by Aesculap.
- Perform sterilization using a dynamic-air-removal (prevacuum) autoclave.
 - Gravity displacement sterilization is not recommended.

- Immediate use (flash) sterilization is not recommended.
- Ensure the sterilizer's maximum load limit is not exceeded when sterilizing multiple sets or devices.
- Do not stack containers as this might prevent the penetration of steam and inhibit drying.
- Refer to the sterilizer manufacturer's instructions and ensure proper installation, calibration, use, and ongoing maintenance.
- The sterilized items should be allowed to cool to room temperature before handling. This allows for safe handling and preventing condensation.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017*.

* Association for the Advancement of Medical Instrumentation (AAMI). Comprehensive guide to steam sterilization and sterility assurance in health care facilities. AAMI ST79:2017. Arlington, VA.

- The following table shows the minimum parameters validated* to achieve a required Sterility Assurance Level (SAL) of 10^{-6} for the system.

Important:

- Sterilization parameters are only valid for devices that have been cleaned per these instructions and are thoroughly dry.
- Sterilization parameters are only valid when the devices are properly housed in the Acumed and Aesculap storage case part numbers identified in the table.
- Please refer to the system IFU for more details on device cleaning, processing, and sterilization methods.

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

Validated Configurations and Parameters:

System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Aculoc 2 Plating System	Case Base: NA Case Lid: NA Trays: 80-0347, 80-0346, 80-0754, 80-0787, 80-0752	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Acutrak 2 Screw System	Case Base: NA Case Lid: NA Trays: AT2-005, AT2-006, AT2-007, AT2-004, 80-1527	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
	Case Lid: 80-0809, 80-0812 Trays: 80-0808, 80-0811	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Acutrak 2 Microextension	Case Base: 80-1527 Case Lid: 80-1534	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Acutrak Fusion	Case Base and Lid: ATF-060	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Anatomic Radial Head System	Trays: TR-0004, TR-0005	Bottom: JN series (perforated Bottom) Lid: JK series	132	25	70
Anatomic Radial Head Solutions	Case Base: NA Case Lid: NA Trays: 80-2002, 80-2003, 80-0365, AT2-006	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Anatomic Radial Head Solutions 2	Case Base: 80-3640/80-3692/80-3693 Case Lid: 80-3641	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Ankle Plating System 3	Case Base: 80-2340 Case Lid: 80-2341	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Arc Wrist Tower	Case Base: NA Case Lid: NA Trays: NA Transfer tray contents into perforated baskets.	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series Aesculap Perforated Basket: JF Series	132	4	60
			134	3	60
Biotrak Resorbable Pin/Nail	Case Base: 80-0485 Case Lid: 80-0486	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Biotrak Resorbable Screw	Case Base: 80-0517 Case Lid: 80-0518	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30

Visit www.aesculap.com/mednet/ifu for the latest version.

System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use Appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Bone Graft System	Case Base: NA Case Lid: NA Transfer tray contents into perforated baskets.	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JF Series	132	4	30
	Case Base: BG-8070 Case Lid: BG-8070	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	134	3	30
6.5/7.3mm Cannulated Screw System	Trays: 80-1825, 80-1827, 80-2020	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
4.0mm Cannulated Screw system	Case Base: 80-3684, 80-3919, 80-3921, 80-3934, 80-4001 Case Lid: 80-3683, 80-3920, 80-3922, 80-3935, 80-4002	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Clavicle Plating System	Trays: 80-0308, 80-0344, 80-0345, 80-0525, 80-0346, 80-0785(optional), 80-0136 (optional), 80-0347 (optional)	Bottom: JN series (perforated Bottom) Lid: JK series	132	25	70
Clavicle Plating System (with updated blue trays)	Trays: 80-0308, 80-0346, 80-0525, 80-3789, 80-3790, 80-3791(optional), 80-0136 (optional), 80-0347 (optional)	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Elbow Plating System	Case Base: NA Case Lid: NA Trays: 80-0675, 80-0676, 80-0677, 80-0678, 80-0629	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Fibula Rod System	Case Base: NA Case Lid: NA Trays: 80-0117, 80-0116	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Forearm Rod System	Trays: ROD-1001 (tray only)	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Fibula and Forearm Nail 2 Base Set	Case Base: 80-2521, 80-3946, 80-3948 Case Lid: 80-2522, 80-3947, 80-3949	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Fibula and Forearm Nail 2 Base Set	Case Base: 80-2521, 80-3946, 80-3948 Case Lid: 80-2522, 80-3947, 80-3949	Bottom: JN series (perforated Bottom) Lid: JK series	134	3	30
Forearm Fracture Solutions	Trays: 80-1915, 80-1916, 80-1922	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Hand Plating System	Case Base: NA Case Lid: NA Trays: 80-0346, 80-0347, 80-1815, 80-1817	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30

System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use Appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Large Acutrak 2 Screw System	Case Base: NA Case Lid: NA Trays: 80-0870, 80-0871, 80-0876, 80-0877	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Lower Extremity Modular System with Ankle and Calcaneal trays OR Forefoot Midfoot Plating Trays	Trays: 80-0427, 80-0428, 80-1501 (Caddy), 80-0744 (Optional), 80-0434 (Optional), 80-0590 (Optional)	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Midshaft Forearm Plate Platter	Trays: 80-0183	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JF Series	132	4	30
Modular Hand System	Case Base and Lid: PL-WF01	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	60
			134	3	60
Polarus 3 Solution	Set 1 Trays: 80-1928, 80-1929, 80-1930 Set 2 Trays: 80-1934, 80-1935, 80-1936	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Polarus Humeral Rod	Trays: MS-9000 (trays only)	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JE Series	132	4	30
Acutrak Screw Removal System	Case Base: 80-0587 Case Lid: 80-0588	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Small Joint Reamer System	Case Base and Lid: MTP-0500	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Small Fragment Base Set	Case Base: 80-2383 Case Lid: 80-2384/ 80-2683	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Total Wrist Fusion Plating System	Case Lid: 80-2123 Trays: 80-0346, 80-1223	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Ulna Shortening Plate (Osteotomy System)	Case Base: 80-0512 Case Lid: 80-0511	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30

OBSOLETE
Visit www.aesculap.com for the latest version.

Symbols Glossary

Symbol	Description
	The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name.

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Sterilization in Aesculap® Rigid Sterilization Containers

EN

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The EN instructions are intended for users in English speaking countries outside the United States and its territories.

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Additional Aesculap Instructions for Use are provided by Aesculap. These instructions can be found by Aesculap Item Number at <http://www.aesculapusa.com/> and must be reviewed prior to use of Aesculap sterile containers. Item Numbers can be found in the table below.

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STERILIZATION

- These products are non-sterile, and must be sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately. Sterilization must be completed in adherence to the additional Instructions for Use provided by Aesculap.
- Perform sterilization using a dynamic-air-removal (prevacuum) autoclave.
 - Gravity displacement sterilization is not recommended.
 - Immediate use (flash) sterilization is not recommended.
- Ensure the sterilizer's maximum load limit is not exceeded when sterilizing multiple sets or devices.

- Do not stack containers as this might prevent the penetration of steam and inhibit drying.
- Refer to the sterilizer manufacturer's instructions and ensure proper installation, calibration, use, and ongoing maintenance.
- The sterilized items should be allowed to cool to room temperature before handling. This allows for safe handling and preventing condensation.
- Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017*.

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System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use Appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Aculoc 2 Plating System	Case Base: NA Case Lid: NA Trays: 80-0347, 80-0346, 80-0754, 80-0787, 80-0752	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Acutrak 2 Screw System	Case Base: NA Case Lid: NA Trays: AT2-005, AT2-006, AT2-007, AT2-004, 80-1527	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
	Case Lid: 80-0809, 80-0812 Trays: 80-0808, 80-0811	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
Acutrak 2 Microextension	Case Base: 80-1527 Case Lid: 80-1534	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Acutrak Fusion	Case Base and Lid: ATF-060	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Anatomic Radial Head System	Trays: TR-0004, TR-0005	Bottom: JN series (perforated Bottom) Lid: JK series	132	25	70
Anatomic Radial Head Solutions	Case Base: NA Case Lid: NA Trays: 80-2002, 80-2003, 80-0365, AT2-006	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Anatomic Radial Head Solutions 2	Case Base: 80-3640/80-3692/80-3693 Case Lid: 80-3641	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Ankle Plating System 3	Case Base: 80-2340 Case Lid: 80-2341	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Arc Wrist Tower	Case Base: NA Case Lid: NA Trays: NA Transfer tray contents into perforated baskets	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series Aesculap Perforated Basket: JF Series	132	4	60
			134	3	60
Biotrak Resorbable Pin/Nail	Case Base: 80-0485 Case Lid: 80-0486	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Biotrak Resorbable Screw	Case Base: 80-0517 Case Lid: 80-0518	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30

Visit www.aesculap.com for the latest version.

System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use Appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Bone Graft System	Case Base: NA Case Lid: NA Transfer tray contents into perforated baskets.	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JF Series	132	4	30
	Case Base: BG-8070 Case Lid: BG-8070	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	134	3	30
6.5/7.3mm Cannulated Screw System	Trays: 80-1825, 80-1827, 80-2020	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
4.0mm Cannulated Screw system	Case Base: 80-3684, 80-3919, 80-3921, 80-3934, 80-4001 Case Lid: 80-3683, 80-3920, 80-3922, 80-3935, 80-4002	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Clavicle Plating System	Trays: 80-0308, 80-0344, 80-0345, 80-0525, 80-0346, 80-0785(optional), 80-0136 (optional), 80-0347 (optional)	Bottom: JN series (perforated Bottom) Lid: JK series	132	25	70
Clavicle Plating System (with updated blue trays)	Trays: 80-0308, 80-0346, 80-0525, 80-3789, 80-3790, 80-3791(optional), 80-0136 (optional), 80-0347 (optional)	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Elbow Plating System	Case Base: NA Case Lid: NA Trays: 80-0675, 80-0676, 80-0677, 80-0678, 80-0629	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Fibula Rod System	Case Base: NA Case Lid: NA Trays: 80-0117, 80-0116	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Forearm Rod System	Trays: ROD-1001 (tray only)	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Forearm Fracture Solutions	Trays: 80-1915, 80-1916, 80-1922	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
Hand Plating System	Case Base: NA Case Lid: NA Trays: 80-0346, 80-0347, 80-1815, 80-1817	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Large Acutrak 2 Screw System	Case Base: NA Case Lid: NA Trays: 80-0870, 80-0871, 80-0876, 80-0877	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Lower Extremity Modular System with Ankle and Calcaneal trays OR Forefoot Midfoot Plating Trays	Trays: 80-0427, 80-0428, 80-1501 (Caddy), 80-0744 (Optional), 80-0434 (Optional), 80-0590 (Optional)	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30

Visit www.aesculap.com for the latest version.

System Name	Case Base and Lid or Tray/Standalone caddy part numbers	Aesculap Part Number (Use Appropriate size containers)	Pre-vacuum Autoclave		
			Exposure Temperature (°C)	Exposure Time (minutes)	Dry Time (minutes)
Midshaft Forearm Plate Platter	Trays: 80-0183	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JF Series	132	4	30
Modular Hand System	Case Base and Lid: PL-WF01	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	60
			134	3	60
Polarus 3 Solution	Set 1 Trays: 80-1928, 80-1929, 80-1930 Set 2 Trays: 80-1934, 80-1935, 80-1936	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Polarus Humeral Rod	Trays: MS-9000 (trays only)	Bottom: JN series (perforated Bottom) Lid: JK series Aesculap Perforated Tray Basket: JE Series	132	4	30
Acutrak Screw Removal System	Case Base: 80-0587 Case Lid: 80-0588	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Small Joint Reamer System	Case Base and Lid: MTP-0500	Bottom: JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Small Fragment Base Set	Case Base: 80-2383 Case Lid: 80-2384/ 80-2683	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Total Wrist Fusion Plating System	Case Lid: 80-2123 Trays: 80-0346, 80-1223	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30
Ulna Shortening Plate (Osteotomy System)	Case Base: 80-0512 Case Lid: 80-0511	Bottom: JK series (Solid Bottom)/ JN series (perforated Bottom) Lid: JK series	132	4	30
			134	3	30

OBSOLETE
 Visit www.aesculapmed.net/ifu for the latest version.

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Fax: +1.503.520.9618

www.acumed.net

These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-76-J | Effective 09-2020 | © 2020 Acumed® LLC

Acumed® is a registered trademark of Acumed, LLC.

Subsidiaries:

Acumed Ltd

Huebner House
The Fairground
Andover
Hampshire UK SP11 0QN
Tel: +44 1264 774450

Acumed Iberica

C/ Alvaro Caballero, 14
28023 Madrid, Spain
Tel: +34 913516357

Acumed Beijing

Room A1206, Horizon International Tower
No. 6, Zhichun Road
Haidian District
100088 Beijing, China
Tel: +86 10 82001303

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Acumed GmbH