



LOGIC™ JR.

PEDIATRIC DISTRACTION SYSTEM



Rethinking Possibilities, Reshaping Lives



LOGIC™ JR.
PEDIATRIC DISTRACTION SYSTEM

- Low profile intra-oral device for use in cases where distraction osteogenesis is indicated
- Utilizes a patented flexible wire that provides access through soft tissue for ease of activation
- Choice of two devices with fixed curvilinear trajectories that mimic the natural growth curve of the mandible

RETHINKING POSSIBILITIES, RESHAPING LIVES

At OsteoMed, we work alongside surgeons who devote their own lives to improving those of others. Our single goal: rethinking the possibilities that have the power to reshape lives.

With a renewed spirit of innovation and a relentless focus on quality, we offer surgeons a wide range of choices for neurosurgical and craniomaxillofacial procedures.




Rethinking Possibilities, Reshaping Lives



UNIQUE • VERSATILE • COMPREHENSIVE



Implant Module




The OsteoMed Logic Jr. Pediatric Mandibular Distraction System is an intraoral bone distractor available in both straight and curvilinear trajectories. The patented curvilinear design allows for the distraction to occur along the natural growth spiral of the mandible. Furthermore, the Logic Jr. Short distractors are a great solution for neonatal mandibular distraction, providing surgeons with a low profile curvilinear design which allows for up to 15mm of distraction length. All other Logic Jr. distractors provide up to 25mm of distraction.




Logic™ Jr. Distractors

52mm Left
(25mm of distraction)




216-1000 LOGIC™ Jr. Distractor,
Left 52mm

52mm Right
(25mm of distraction)




216-1001 LOGIC™ Jr. Distractor,
Right 52mm

52mm Straight Left/Right
(25mm of distraction)




216-1002 LOGIC™ Jr. Distractor,
Straight (Left/Right)

52mm Left Short
(15mm of distraction)



216-1003 LOGIC™ Jr. 52mm Left
Short Distractor

52mm Right Short
(15mm of distraction)



216-1004 LOGIC™ Jr. 52mm Right
Short Distractor

1.2mm Screws

1.2mm Standard Screws



204-1202 1.2mm x 2mm Screw
204-1203 1.2mm x 3mm Screw
204-1204 1.2mm x 4mm Screw
204-1206 1.2mm x 6mm Screw
204-1208 1.2mm x 8mm Screw
204-1210 1.2mm x 10mm Screw
204-1212 1.2mm x 12mm Screw

1.2mm AutoDrive® Screws



211-1203 1.2mm x 3mm Auto-Drive® Screw
211-1204 1.2mm x 4mm Auto-Drive® Screw
211-1205 1.2mm x 5mm Auto-Drive® Screw

1.5mm Safety Screws



204-1503 1.5mm x 3mm
Safety Screw
204-1506 1.5mm x 6mm
Safety Screw

Activation Wires



216-0302 Activation Wire, 62mm



216-0306 Activation Wire, 72mm



216-0307 Activation Wire, 82mm



216-0305 Silicone Tubing

Instruments



220-0055 Short Pointed Cannula



220-0019 TAPERLOCK™ Screwdriver Body



220-0027 Small Grasping Forceps



216-0102 Distraction Tool



220-0049 Plate Bending Forceps



220-0028 Plate Cutter



220-0056 Trocar, Short Pointed Cannula



220-0058 Cheek Retractor



220-0140 Drill Guide, Neutral, Short Pointed Cannula



216-0227 2mm Spacer



216-0228 4mm Spacer



220-0085 1.2mm AUTO-DRIVE® Screw VISIDISK™



220-0071 1.2mm Screw VISIDISK™



Rethinking Possibilities, Reshaping Lives



OSTEOMED

3885 Arapaho Rd.
Addison, TX 75001

Customer Service: 800.456.7779
Customer Service Fax: 800.390.2620

www.osteomed.com

P/N 030-1439 Rev.F