




ExtremiLock™ Foot Plating System

Lisfranc Plating System

Offering **five plate styles and 25 left and right-specific plate options** in various sizes, providing an ideal fit for diverse patient anatomies



Lisfranc Targeting Guide
improves accuracy for
Lisfranc screw placement




Fix or Fuse? A New Comprehensive Solution for Both Approaches, Customizing Patient Care

Lisfranc Plating System

Designed in conjunction with Lawrence Fallat, DPM, J. Robert Faux, MD, and Nirmal Tejawani, MD, the Lisfranc Plating System provides fixation for Lisfranc fractures-dislocations and fusions of the tarsal-metatarsal joints. This system offers the most complete selection of plates for a dorsal, dorsomedial, and medial approach. The plate design allows for compression and stabilization of the TMT joint while allowing surgeons to visualize the affected area during the healing process.


The system offers a family of five plates that include 25 left and right-specific plate options in several sizes to fit various patient anatomies. The plates are low-profile and are anatomically contoured to fit the TMT joint.



System Features

Lisfranc Targeting Guide

Facilitates accurate Lisfranc screw fixation placement and trajectory for intermetatarsal instability cases





2nd—3rd Dorsal Plate



Dorsomedial Plate



View More X-rays

go.acumed.net/LF-Xrays

Universal Screw Holes

All threaded holes accept 2.7 mm, 3.5 mm, and/or 4.0 mm locking and nonlocking screws in any hole, and up to 20° of variable-angle locking in any direction (40° conical)



Medial Fixation

Distal hole trajectory allows use of interfragmentary screws and/or flexible fixation implant for intermetatarsal instability fixation



Anatomically Contoured

Multiple Distal Fixation Points

Optimizing load distribution, screw purchase, and minimizing stress fractures

Contact Your Local Sales Rep for More Information



1st—2nd Dorsal Plates



2nd—3rd Dorsal Plates



Dorsomedial Plates



Stabilization Plates



Universal Medial Plate

Learn More About ExtremiLock™ Foot Plating System

The Most Complete Selection of Lisfranc
Fractures/Dislocations and Fusions

go.acumed.net/Lisfranc



go.acumed.net/Elock-FPS

OsteoMed ExtremiLock Foot
Plating System



go.acumed.net/Cann-Mini-Small

OsteoMed ExtremiFix Cannulated
Screw System—Mini & Small



www.acumed.net

Acumed USA Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
+1.888.627.9957

OsteoMed USA Campus
3885 Arapaho Road
Addison, TX 75001
+1.800.456.7779

Acumed Iberica Campus
C. Proción, 1
Edificio Oficor
28023 Madrid, Spain

FNA10-04-A | Effective: 2023/09 | © 2023 Acumed® LLC

OsteoMed and ExtremiLock™ are registered trademarks of OsteoMed LLC.
Acumed® is a registered trademark of Acumed LLC.